

# Tuttnauer

Your Sterilization & Infection Control Partners

## Sterilization Without Compromise

- Expertise – Specializing in Sterilization
- Experience - Over 85 Years
- Choice - Largest Line to Choose from

**ALL YOU NEED TO DO IS CHOOSE**



Fully Automatic Autoclaves



Elara11 - Class B  
Pre and Post Vacuum Autoclaves



Large Capacity Autoclaves



Manual Autoclaves



Compact Autoclaves



Water Distillers



Ultrasonic Cleaners and Tablets



Autoclave Cleaner

Model	Chamber Dimensions Ø x L in(mm)	Volume	Overall Dimensions DxWxH in(mm)	Standard Cassette Capacity	Tray Dimensions DxWxH in(mm)	No. of Trays	Standard Unwrapped Cycle Time	Voltage(V) Freq.(Hz)	Power(W) Current(A)	Shipping Weight lbs(kg)
-------	---------------------------------------	--------	---------------------------------------	----------------------------	------------------------------------	--------------	-------------------------------	-------------------------	------------------------	----------------------------

### Fully Automatic Autoclaves

<b>EZ9</b>	9" x 18" (230 x 470)	5gal/19L	21.5"x20"x14.4" (545x510x365)	2 full and 2 half	16.3"x 6.7"x0.8" (415x170x20)	3	Cold: 27 mins Hot: 13 mins	120V 50/60	1400W 12A	95 (43)
<b>EZ10</b>	10" x 19" (254 x 475)	6gal/23L	21.5"x20"x14.4" (545x510x365)	3 full and 3 half	16.3"x 6.7"x0.8" (415x170x20)	4	Cold: 30 mins Hot: 14 mins	120V 50/60	1400W 12A	95 (43)
<b>EZ10k</b>	10" x 19" (254 x 475)	6gal/23L	21.5"x20"x14.4" (545x510x365)	3 full and 3 half	16.3"x 6.7"x0.8" (415x170x20)	4	Cold: 21 mins Hot: 11 mins	230V* 50/60	2200W 10A	95 (43)
<b>3870EA</b>	15" x 30" (384 x 760)	22gal/85L	34.5"x26"x20.6" (875x665x525)	15 full	26.6"x11.3"x1" (675x286x25) 26.6"x13.8"x1" (675x350x25)	2	Cold: 32 mins Hot: 23 mins	230V* 50/60	3000W 13.6A	213 (97)

### Manual Autoclaves

<b>Valueklave 1730</b>	7" x 13" (170 x 340)	2gal/7.5L	17.9"x17.4"x12" (455x440x305)	2 half	11.6"x4.7"x0.8" (295x120x20)	3	Cold: 16 mins Hot: 11 mins	120V 50/60	1350W 11A	50 (23)
<b>2340M</b>	9" x 18" (230 x 470)	5gal/19L	21.5"x20"x14.4" (545x510x365)	2 full and 2 half	16.3"x6.7"x0.8" (415x170x20)	3	Cold: 27 mins Hot: 13 mins	120V 50/60	1400W 12A	85 (39)
<b>2540M</b>	10" x 19" (254 x 475)	6gal/23L	21.5"x20"x14.4" (545x510x365)	3 full and 3 half	16.3"x6.7"x0.8" (415x170x20)	4	Cold: 30 mins Hot: 14 mins	120V 50/60	1400W 12A	85 (39)
<b>2540MK</b>	10" x 19" (254 x 475)	6gal/23L	21.5"x20"x14.4" (545x510x365)	3 full and 3 half	16.3"x6.7"x0.8" (415x170x20)	4	Cold: 21 mins Hot: 11 mins	230V* 50/60	2200W 10A	85 (39)
<b>3870M</b>	15" x 30" (384 x 760)	22gal/85L	34.5"x26"x20.6" (875x665x525)	15 full	26.6"x11.3"x1" (675x286x25) 26.6"x13.8"x1" (675x350x25)	2	Cold: 32 mins Hot: 23 mins	230V* 50/60	3000W 13.6A	209 (95)

### Pre and Post Vacuum (Class B)

<b>Elara11</b>	11" x 19.8" (280 x 504)	7.5 gal/28.5L	25.4"x20.8"x17.4" (645x530x444)	4 full and 4 half	16.3"x6.7"x0.8" (415x170x20)	5	Cold: 15 mins Hot: 15 mins	230V* 50/60	2300W 10A	184 (84)
----------------	----------------------------	---------------	------------------------------------	----------------------	---------------------------------	---	-------------------------------	----------------	--------------	-------------

\*Requires a voltage between 220V and 235V for proper operation. When applicable a buck/boost transformer (0.5KVA) may be required.

[www.tuttnauerusa.com](http://www.tuttnauerusa.com)

Tuttnauer USA Co. Ltd., 25 Power Drive, Hauppauge, NY 11788  
Tel: (800) 624 5836, (631) 737 4850, Fax: (631) 737 0720  
Email: [info@tuttnauerUSA.com](mailto:info@tuttnauerUSA.com), [www.tuttnauerUSA.com](http://www.tuttnauerUSA.com)

